

## Nucleic Acid Labeling and Delivery Tools for Neurobiology

Aaron G. Loomis, Jennifer L. Duzeski, Laura K. Juckem, Kira J. Machnik, Shyla L. Tessmer, Mary-Anne V.

Watt, James E. Hagstrom. Mirus Corporation, 2002.

### Introduction

To address the needs of neuroscientists, Mirus has developed a variety of optimized reagents for labeling and transfection of nucleic acids. Mirus' **TransIT<sup>®</sup>-Neural™ Reagent** is a plasmid transfection reagent specifically designed to achieve the highest transfection efficiency, while maintaining cell viability, in a wide range of neural cell lines. Mirus' popular **TransIT-TKO<sup>®</sup> Transfection Reagent** efficiently delivers siRNA, enabling knock out of target gene expression in a wide variety of mammalian cell types, including those that are of neural origin. Mirus also provides **Label IT<sup>®</sup> Tracker™ Kits** and **Label IT<sup>®</sup> siRNA Tracker™ Kits** for one-step direct covalent labeling of nucleic acids in a non-destructive manner, allowing for intracellular tracking studies, where the original function of the nucleic acid is generally maintained. The **Label IT<sup>®</sup> Tracker™ Kits** and **Label IT<sup>®</sup> siRNA Tracker™ Kits** are available in a variety of labels (Cy™ 3, Cy™ 5, CX-Rhodamine, TM-Rhodamine, Fluorescein, and Biotin), which are ideal for tracking nucleic acids with one or more fluorophores. In this report, we demonstrate how these reagents, when used individually or in concert, facilitate effective neurobiological studies in gene delivery and regulation.

### Methods

#### Preparation of Labeled pDNA

Plasmid DNA was labeled with the indicated **Label IT<sup>®</sup> Tracker™ Reagent** at a ratio of 1 µl reagent per 1 µg of DNA and incubated for 1 hour at 37°C, as described in the **Label IT<sup>®</sup> Tracker™** protocol. The sample was purified by ethanol precipitation, and the labeled DNA pellet was resuspended in sterile, molecular-biology grade water and

quantified. Labeled DNA was stored at -20°C, protected from light.

#### Transfection of pDNA

Transfection complexes were prepared by adding 1-16 µl (depending on cell type) of the **TransIT<sup>®</sup>-Neural™ Transfection Reagent** per 1 µg of plasmid DNA to 100 µl of serum-free media, vortexing and incubating at room temperature for 5-10 minutes. A firefly luciferase plasmid (1 µg per well of a 12-well plate) was added to the diluted **TransIT<sup>®</sup>-Neural™ Reagent**, mixed by gentle pipetting and incubated at room temperature for an additional 5-10 minutes. Complexes were added dropwise to 50-70% confluent cells in their complete media and mixed by gentle rocking. Following the addition of the complexes, the cells were incubated at 37°C for 24-48 hours and assayed for gene expression.

#### Preparation of siRNA Duplexes

An RNA oligonucleotide directed against the coding region for firefly luciferase and its reverse complement were synthesized and PAGE purified. Anti-firefly luciferase siRNA duplexes were formed by hybridizing the reverse complementary RNA oligonucleotides in 100 mM NaCl, 50 mM Tris pH 8 buffer. These were then diluted in the same buffer to a working concentration of 1 µM.

#### Preparation of Labeled siRNA

Duplex siRNA was labeled with the indicated **Label IT<sup>®</sup> siRNA Tracker™ Reagent** at a ratio of 1 µl reagent per 1 µg of siRNA and incubated for 1 hour at 37°C, as described in the **Label IT<sup>®</sup> siRNA Tracker™** protocol. The sample was purified by ethanol precipitation, and the labeled siRNA pellet was resuspended in the siRNA Dilution Buffer (provided in the kit) and quantified. Labeled siRNA was stored at -20°C, protected from light.

### Transfection of siRNA Duplexes

Transfection complexes were prepared by adding 1-8  $\mu$ l of *TransIT*-TKO<sup>®</sup> Reagent to 50-200  $\mu$ l (volumes depend on plate size) of serum-free media, vortexing and incubating the mixture at room temperature for 5-10 minutes. Anti-firefly luciferase siRNA duplex (5-50 nM final concentration in the well) was added to the mixture, mixed by gentle pipetting and incubated at room temperature for an additional 5-10 minutes. Complexes were added dropwise to 50-70% confluent cells in their complete media and mixed by gentle rocking. Cells were incubated at 37°C for 24-48 hours and assayed for gene expression.

### Simultaneous Transfection of pDNA and siRNA Duplexes


The transfection complexes were prepared in the specified order. One to two  $\mu$ l of *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent were added to 50  $\mu$ l of serum-free media, vortexed and incubated at room temperature for 5-10 minutes. Two plasmid DNAs, one encoding sea pansy luciferase (5 ng) and the other firefly luciferase (500 ng) were added to the diluted *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent by gentle pipetting and incubated at room temperature for 5-10 minutes. To the same mixture, 0.5-2  $\mu$ l of *TransIT*-TKO<sup>®</sup> Reagent was added and mixed by gentle pipetting and incubated at room temperature for 5-10 minutes. Anti-firefly siRNA duplex (5-10nM) was then added and mixed by gentle pipetting and incubated at room temperature for 5-10 minutes. This entire mixture was added dropwise to 50-70% confluent cells (cultured in 24-well plates) in their complete media and mixed by gentle rocking. Cells were incubated at 37°C for 24 hours and assayed for gene inhibition.

## **Results and Discussion**

### High Efficiency Neural Cell Transfection

*TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent has been optimized to successfully deliver plasmid DNA to many neural cell lines, including C6, Daoy, DBTRG-05MG, DI-TNC1, HCN-1A, Neuro-2a, PC-12, SK-N-MC, SVG p12 and human astrocytes. Neuro-2a cells, a common neuroblastoma cell line in neurobiology research, was used to demonstrate the ability of *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent to deliver EGFP plasmid DNA. Forty-eight hours post-transfection, the cells were fixed and mounted for fluorescent microscopy (Figure 1). At least

75% of the Neuro-2a cells were successfully transfected, while cell viability was maintained.




**Figure 1. Expression of pEGFP in Neuro-2a Cells**

Neuro-2a cells transfected with 4  $\mu$ l *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent and 2  $\mu$ g pEGFP (Clontech) per well of a 6-well plate and incubated in complete media. The image on the left was acquired 48 hours post-transfection on a Zeiss Axiocam. The image on the right was acquired 48 hours post-transfection on a Zeiss LSM 510 Confocal Microscope.

### Plasmid DNA Tracking and Expression in Neural Cells

Labeled plasmid DNA can be introduced into neural cells in culture by a variety of methods. Transfections using *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent are characterized by low toxicity and high efficiency. As seen in Figure 2, *Label IT*<sup>®</sup> Tracker<sup>™</sup> Cy<sup>™</sup>5 labeled EGFP plasmid DNA can be seen in all cells 24 hours post-transfection; EGFP expression can be visualized simultaneously. See Mirus' *Label IT*<sup>®</sup> Tracker<sup>™</sup> Intracellular Nucleic Acid Localization Kits Technical Report for more details on tracking and expression of plasmid DNA.




**Figure 2. Tracking and Expression of Plasmid DNA in DI-TNC1 Cells**

DI-TNC1 astrocytes were transfected with 1  $\mu$ g *Label IT*<sup>®</sup> Tracker<sup>™</sup> Cy<sup>™</sup>5 labeled pEGFP and 3  $\mu$ l *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent per well of a 6-well plate in complete media. The blue represents the Cy<sup>™</sup>5 fluorescence, and the green represents pEGFP expression. Images were acquired at 24 hours post-transfection using a Zeiss LSM 510 Confocal Microscope.

### Knock Out of Target Gene Expression

siRNA complexed with *TransIT*-TKO<sup>®</sup> Transfection Reagent has the ability to selectively knock out target gene

expression without affecting the expression of a similar non-targeted gene. To illustrate this, firefly luciferase and sea pansy luciferase expression vectors were co-transfected into a variety of cell types with *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent. To the same tube, *TransIT*-TKO<sup>®</sup> Reagent/anti-firefly luciferase siRNA complexes were added. Using 5-10 nM siRNA, greater than 85% knock out of firefly luciferase expression was observed in all cell lines tested at 24 hours when compared to the expression of the untargeted sea pansy luciferase (Figure 3).




**Figure 3. Target Gene Expression Knock Out in Neural Cell Lines**

Efficiency of knock out in various neural lines using the *TransIT*-TKO<sup>®</sup> Transfection Reagent with 5-10 nM anti-firefly luciferase siRNA duplexes (Dharmacon Research, Inc). Cell lines were first transfected with plasmids expressing firefly and sea pansy luciferase using the *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent. Twenty-four hours post-transfection, cell lysates were assayed for both firefly and sea pansy luciferase expression using Promega's Dual-Luciferase<sup>®</sup> Reporter Assay System.


#### Tracking of siRNA Duplexes

In order to visualize the subcellular location of siRNA after delivery by *TransIT*-TKO<sup>®</sup> Reagent, *Label IT*<sup>®</sup> siRNA Tracker<sup>™</sup> Cy<sup>™</sup>3 labeled siRNA duplexes were transfected into human astrocytes (Figure 4) and C6 glioma cells (Figure 5) and then visualized using confocal microscopy. The majority of the siRNA is cytoplasmic. Both diffuse and punctate signals were observed, indicating that there are likely at least two populations of delivered labeled siRNA molecules; one which is free in the cytoplasm itself and the other which is contained within endosomes.



**Figure 4. siRNA Tracking in Human Astrocytes**

Human astrocytes (ScienCell Research Laboratories) were seeded on poly-D-lysine coated coverslips in a 24-well plate and incubated overnight. One  $\mu$ l of *TransIT*-TKO<sup>®</sup> Transfection Reagent was diluted in 50  $\mu$ l of serum-free media and 50 nM (final concentration in well) of *Label IT*<sup>®</sup> siRNA Tracker<sup>™</sup> Cy<sup>™</sup>3 labeled siRNA duplex was added to the diluted *TransIT*-TKO<sup>®</sup> Transfection Reagent for complex formation. Subsequently, these complexes were added to 30% confluent human astrocytes in their complete growth media. Confocal image of the single astrocyte was acquired 24 hours post-transfection on a Zeiss LSM 510 Confocal Microscope.



**Figure 5. siRNA Tracking in C6 Cells**

C6 glioma cells were seeded on poly-D-lysine coated coverslips in a 6-well plate and incubated overnight. Eight  $\mu$ l of *TransIT*-TKO<sup>®</sup> Transfection Reagent were diluted in 200  $\mu$ l of serum-free media and 50 nM (final concentration in well) of *Label IT*<sup>®</sup> siRNA Tracker<sup>™</sup> Cy<sup>™</sup>3 labeled siRNA duplex was added to the diluted *TransIT*-TKO<sup>®</sup> Transfection Reagent for complex formation. Subsequently, these complexes were added to 45% confluent C6 cells in their complete growth media. Confocal image was acquired 20 hours post-transfection on a Zeiss LSM 510 Confocal Microscope.

## Summary

Mirus' *TransIT*<sup>®</sup>-Neural<sup>™</sup> and *TransIT*-TKO<sup>®</sup> Transfection Reagents allow for highly efficient transfection of plasmid DNA and siRNA, respectively, in a variety of neural cell lines. These reagents include optimized protocols which result in high cellular viability and are most effective when carried out in serum containing media, with no media changes or additions. The *Label IT*<sup>®</sup> Tracker<sup>™</sup> and *Label IT*<sup>®</sup> siRNA Tracker<sup>™</sup> Kits allow for the direct covalent labeling of plasmid DNA and siRNA. These tools give researchers the opportunity to assess transfection efficiencies on their particular neural cell line, and perform intracellular localization studies, where the original function of the nucleic acid is generally maintained. The combination of these products provides neurobiologists with simple methods to reproducibly study gene transfer and expression as well as target gene knock out in neural cell lines.

## Product Protocols (available at [www.genetransfer.com](http://www.genetransfer.com))

- *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent Protocol #ML022
- *TransIT*-TKO<sup>®</sup> Transfection Reagent Protocol #ML015
- *Label IT*<sup>®</sup> Tracker<sup>™</sup> Intracellular Nucleic Acid Localization Kits Protocol #ML018
- *Label IT*<sup>®</sup> siRNA Tracker<sup>™</sup> Intracellular Localization Kits Protocol #ML023

## Technical Reports (available at [www.genetransfer.com](http://www.genetransfer.com))

- *TransIT*-TKO<sup>®</sup> Short Interfering RNA (siRNA) Transfection Reagent (Technical Report #MTR001)
- Intracellular Localization and Expression of Labeled Plasmid DNA using *Label IT*<sup>®</sup> Tracker<sup>™</sup> Nucleic Acid Labeling Kits (Technical Report #MTR002)

Mirus Product References are available at [www.genetransfer.com](http://www.genetransfer.com).

For Customer and Technical Support  
contact Mirus at:  
**888.530.0801 or 608.441.2852**  
[www.genetransfer.com](http://www.genetransfer.com)

## Ordering Information

### *TransIT*<sup>®</sup>-Neural<sup>™</sup> Transfection Reagent

<u>Product Name</u>	<u>Quantity</u>	<u>Cat. No.</u>
<i>TransIT</i> <sup>®</sup> -Neural <sup>™</sup> Transfection Reagent	1 ml	MIR 2140
<i>TransIT</i> <sup>®</sup> -Neural <sup>™</sup> Transfection Reagent	0.4 ml	MIR 2144
<i>TransIT</i> <sup>®</sup> -Neural <sup>™</sup> Transfection Reagent	5 x 1 ml	MIR 2145
<i>TransIT</i> <sup>®</sup> -Neural <sup>™</sup> Transfection Reagent	10 x 1 ml	MIR 2146

### *TransIT*-TKO<sup>®</sup> Transfection Reagent

<u>Product Name</u>	<u>Quantity</u>	<u>Cat. No.</u>
<i>TransIT</i> -TKO <sup>®</sup> Transfection Reagent	1 ml	MIR 2150
<i>TransIT</i> -TKO <sup>®</sup> Transfection Reagent	0.4 ml	MIR 2154
<i>TransIT</i> -TKO <sup>®</sup> Transfection Reagent	5 x 1 ml	MIR 2155
<i>TransIT</i> -TKO <sup>®</sup> Transfection Reagent	10 x 1 ml	MIR 2156

### *Label IT*<sup>®</sup> Tracker<sup>™</sup> Intracellular Nucleic Acid Localization Kits

<u>Product Name</u>	<u>Quantity</u>	<u>Cat. No.</u>
<i>Label IT</i> <sup>®</sup> Tracker <sup>™</sup> Cy <sup>™</sup> 3	*	MIR 7010
<i>Label IT</i> <sup>®</sup> Tracker <sup>™</sup> Cy <sup>™</sup> 5	*	MIR 7011
<i>Label IT</i> <sup>®</sup> Tracker <sup>™</sup> CX-Rhodamine	*	MIR 7012
<i>Label IT</i> <sup>®</sup> Tracker <sup>™</sup> TM-Rhodamine	*	MIR 7013
<i>Label IT</i> <sup>®</sup> Tracker <sup>™</sup> Fluorescein	*	MIR 7015
<i>Label IT</i> <sup>®</sup> Tracker <sup>™</sup> Biotin	*	MIR 7014

\*Each kit contains *Label IT*<sup>®</sup> Tracker<sup>™</sup> Reagent, Tracker Reconstitution Solution, 10X Labeling Buffer A, and *TransIT*<sup>®</sup>-LT1 Transfection Reagent. Each kit provides sufficient reagents to label 50 to 200 µg of plasmid DNA and perform at least 50 transfections

### *Label IT*<sup>®</sup> siRNA Tracker<sup>™</sup> Intracellular Localization Kits

<u>Product Name</u>	<u>Quantity</u>	<u>Cat. No.</u>
<i>Label IT</i> <sup>®</sup> siRNA Tracker <sup>™</sup> Cy <sup>™</sup> 3	**	MIR 7200
<i>Label IT</i> <sup>®</sup> siRNA Tracker <sup>™</sup> Cy <sup>™</sup> 5	**	MIR 7201
<i>Label IT</i> <sup>®</sup> siRNA Tracker <sup>™</sup> CX-Rhodamine	**	MIR 7202
<i>Label IT</i> <sup>®</sup> siRNA Tracker <sup>™</sup> TM-Rhodamine	**	MIR 7203
<i>Label IT</i> <sup>®</sup> siRNA Tracker <sup>™</sup> Fluorescein	**	MIR 7205
<i>Label IT</i> <sup>®</sup> siRNA Tracker <sup>™</sup> Biotin	**	MIR 7204

\*\*Each kit contains *Label IT*<sup>®</sup> siRNA Tracker<sup>™</sup> Reagent, *Label IT*<sup>®</sup> Reconstitution Solution, 10X Labeling Buffer A, siRNA Dilution Buffer, and *TransIT*-TKO<sup>®</sup> Transfection Reagent. Each kit provides sufficient reagents to label 50 µg of siRNA and perform up to 500 transfections in 24-well plates.

*TransIT*, *TransIT*-TKO, and *Label IT* are registered trademarks of Mirus Corporation. *TransIT*-Neural, Tracker, and Gene Transfer Specialists are trademarks of Mirus Corporation. Mirus Transfection Reagents are covered by United States Patent No. 5,744,335; 5,965,434; 6,180,784; 6,383,811 and patents pending. *Label IT* Reagents are covered by U.S. Patent No. 6,262,252.

Cy3 and Cy5 are trademarks of Amersham Biosciences. Cy3 and Cy5 products or portions thereof are manufactured under license from Carnegie Mellon University and are covered by U.S. Patent No. 5,268,486.

© 2002, Mirus Corporation. All rights reserved.